


## elmex® interX toothbrushes

## Study result

## Cleaning performance of elmex® interX toothbrushes

The cleaning performance of the elmex® interX toothbrushes was tested in a standardised in vitro model. A conventional flat-profile toothbrush served as control brush.



The elmex® interX toothbrushes show an excellent interproximal access efficacy, allowing more plaque removal in the difficult to reach interdental areas. Short head toothbrushes are even more effective.

## Reference

Yankell SL, Shi X: Laboratory evaluations of elmex® interX toothbrushes for interproximal access efficacy and posterior tooth surface cleaning. J Clin Dent 13 (2002), 249-252



CH-4142 Münchenstein, Switzerland www.gaba.com